

**ClariVein OC- Procedure List**

- ClariVein unit- 45cm, 65cm or 85cm
- The bed needs to be able to place the patient in the supine and prone position
- Sterile table cover
- Sterile probe cover
- Sterile Ultrasound gel
- U-drape. Or appropriate drape for patient (193 cm x 305 cm Split 10 cm x 102 cm)
- Prep to sterile the leg
- Local anaesthetic
- NEEDLE, 23G X 30 MM
- NEEDLE, BLUNT FILTER 18G X 40 MM
- MAK 4Fr or 5Fr- echo enhanced needle
- Sharps holder
- GAUZE 10 X 10CM 12PLY
- MEPORE, FILM AND PAD, 5CM X 7CM
- 2 gallipots- 1 for saline, 1 for mixed Fibrovein solution
- 2 syringes 10ML
- Compression Stocking

**Sclerosant- Fibrovein or Polidocanol****Fibrovein**

- 1% and 3%
- Saline
- 3% Fibrovein is diluted with saline to make 1.5%

**Polidocanol**

- 2% Polidocanol

The amount of Fibrovein or Polidocanol depends on the length and average diameter of the saphenous vein.

See dosage chart for reference

